


Sittan

Best Choice for Bed Side Rehabilitation



The sitting position Without Back Support



PARAMOUNT BED

What is the sitting position "Without Back Support" ?

Bad effect from the bedridden

- The risk of pneumonia because of unsatisfied pulmonary ventilation
- The risk of postural hypotension or DVT from lowering body fluid and blood stagnation
- Muscle weakness by joint stiffness or loss of exercise

Start the sitting position "Without Back Support" = First stage of rehabilitation as early as possible

Triple necessary keys to form the sitting position "Without Back Support"

1



Upright sitting without back support

2




Keep the neck by oneself

3



Keep the feet flat on the floor



Note: We recommend to use the Sittan with low height bed to enhance the effectiveness of the position.

Specifications at a glance



Features

- Height adjustable to fit the patient's body
- Length adjustable for body size
- Posture support belt
- Arm support cushion
- Foldable structure for easy storing

Product code		KF-890
Description		Rehabilitation Table
Dimensions	Width x Length	680mm x 770mm
	Height of table top	550mm to 800mm
	Length adjustable range	770mm to 950mm
	Angle of table top	0 - 90°
Material	Caster	Front: 50mm Back: 30mm
	Body	Steel
	Table top	Melamine faced particle board



PARAMOUNT BED

PARAMOUNT BED CO., LTD.
 2-14-5 Higashisuna Koto-Ku, Tokyo 136-8670, Japan
 TEL : +81 (0)3 3648 2961 FAX: +81 (0)3 3648 2420
<http://www.paramount.co.jp/english/>