

Aerogen®



- Refill medication without breaking the circuit
- 4 times more medication delivered to patients lungs as aerosol size is (1-5 μm)
- Can be used with all medications for inhalation
- No Heating or degradation of medication.
- Minimal Residual Volume (<0.1ml for 3ml Dose)
- Operates with ventilated and non-ventilated patients when coupled with the Aerogen Ultra.



- can be placed at the Wye or at the humidifier
- Single patient use for upto 28 days intermittent or 7 days continuous.(Solo only)
- Multiple patient use. Autoclavable at 132C-135C (Pro only)
- Silent operation

