



Aeon7200A

Basic Configuration

Anaesthesia Machine

CE 0123

AEOMED
Reliable Quality Thoughtful Service

- Large LCD screen, various monitoring waveform and parameters in different colors
- IPPV, Pressure mode, SIMV modes can easily meet the ventilation requirement
- Integrated CO₂ absorption system equipped with NO WATER VAPOR CONDENSATION sensor and fast-release APL valve
- 4-tube of flow meters for 2-gas sources meet long operation requirement
- Optional AGSS provides more protection
- Optional multi-gas monitoring
- With CE certificate, meets EU clinical safety requirement

Technical Specification

Physical & Environmental Specification		• Setting:	Navigator knob and touch screen
• Dimensions(H×W×D):	1356×845×620mm	• Tidal volume:	20~1500ml
• Weight:	110kg (basic unit)	• Ventilation frequency:	4~100bpm
• Storage temperature:	-20~+55°C	• I:E ratio:	4:1~1:8
• Working temperature:	+10~+40°C	• PEEP:	0~30 cmH ₂ O (optional)
• Storage humidity:	≤93%	• Inspiratory plateau:	OFF, 5%~60%
Electrical Specification		• Flow trigger:	2L/min~30L/min
• Mains:	100~240V, 50/60Hz	• Pressure trigger:	0~20 cm H ₂ O
• Battery:	24V DC, minimum 120 min	• Pressure range:	5~70 cm H ₂ O
• Mains outlet:	3	• Monitoring display:	Display: 8.4" color touch screen
Pneumatic Specification			Pressure: P _{peak} , P _{plat} , P _{mean} , PEEP
• Gas supply:	0.28~0.6MPa		Volumes: MV, Vt
• Flowmeters:	O ₂ : 0~1L/min; 1~10L/min		Breath rate
	N ₂ O: 0~1L/min; 1~10L/min		Compliance
• Gas system:	Oxygen supply failure alarm		Gas: O ₂ , CO ₂ , Agent (all optional)
	Oxygen supply failure protection		Graph display: P-t, F-t waveforms
	O ₂ flush: 25~75L/min	• Alarm:	High Airway Pressure/Low Airway Pressure
• Working mode:	Closed, Semi-closed, Semi-open		High Minute Volume/Low Minute Volume
• Driven mode:	Pneumatically driven and electronically controlled		High FiO ₂ /Low FiO ₂ (when FiO ₂ function is enabled)
• Safety valve:	≤10kPa		Power failure
• Operating mode:	Man. / Vent.		O ₂ Supply failure, Apnea,
• Cylinder yokes:	Optional		Alarm Silence: ≤120 seconds
Ventilator Specification		Vaporizer	
• Patient type:	Adult, pediatric	• Support 2 vaporizers (Selectatec® with interlock, optional standby vaporizer parking holder)	
• Ventilation modes: IPPV, PLV, SIMV, Spont (with apnea backup), Manual, Standby		• Agent type: Halothane, Enflurane, Isoflurane, Sevoflurane	
		• Filling type: Pour-fill, Key-fill, Quik-fil®	

Remark: Above configurations include standard and option. Please check price with your Aeonmed sales representative.

AEONMED
Reliable Quality Thoughtful Service

HQ: No.4 Hangfeng Road, Fengtai Science Park, Fengtai District, Beijing, China (100070)
Science Park: No. 10, Chaobai Street, Yanjiao Development Zone, Sanhe City, Hebei Province, China (065201)
TEL: +86-10-5841 1198 FAX: +86-10-6371 8989
Http: //www.aeonmed.com E-mail: marketing@aeonmed.com

Information contained in this leaflet is correct at the date of publication. The policy of Beijing Aeonmed Co., Ltd. is one of the continued improvement to its products. Because of this policy Beijing Aeonmed Co., Ltd. reserves the right to make any change, which may affect the information in this leaflet without giving prior notices.

AAM7200AL-1703